


Genotyping *Alk1^{tm2.2Spo}*

(*Acvr1^{tm2.2Spo}*; *Alk1^{1f}* or *Alk1^{1loxP}* = KO; *Alk1^{2f}* = floxed, no neo; *Alk1^{2f}* = floxed, with neo)

MGI: 4398902

Reference: Park SO, et al., *ALK5- and TGFBR2-independent role of ALK1 in the pathogenesis of hereditary hemorrhagic telangiectasia type 2. Blood.* 2008 Jan 15;111(2):633-42



Schematic diagram of the *Alk1* wild-type allele, *Alk1*-conditional targeting vector, and *Alk13loxP*, *Alk12loxP*, and *Alk11loxP* alleles. Exons and loxP sequences are indicated by boxes and arrowheads, respectively. Locations of primer pairs used for amplifying specific regions containing a loxP sequence are also indicated. We generated a conditional *Alk1* knockout allele in which 3 loxP sequences were inserted into intronic regions of the *Alk1* gene. The 5 most loxP sequence and the neomycin cassette flanked by loxP sequences were inserted into the third and the sixth introns, respectively; thus, exons 4-6 could be deleted by Cre-mediated recombinations (Figure 2A). Because exon 5 encodes the transmembrane (TM) domain, the deletion of exons 4-6 would produce a null allele.

PRIMERS:

Alk1 (2f): - conditional allele

oJDW 703 (m*Alk1*-ex5FWD): 5'-cctggacagcgactgtactac

oJDW 704 (m*Alk1*-in6REV): 5'-gccccattgctctcctcaaac

WT = 350 bp

FLOX = 400 bp

Alk1 (1f) -after Cre mediated deletion (aka Null or KO allele)

oJDW 705 (m*Alk1*-in3FWD4): 5'-CAGCACCTACATCTTGGGTGGAGA

oJDW 706 (m*Alk1*-in6REV, same as oJDW 704!): 5'-GCCCCATTGCTCTCCTCAAAC

KO=390bp

REACTION RECIPE:

9.375 ul ddH₂O

1.25 ul 10x Buffer

0.25 ul 10mM dNTPs

0.25 ul JDW 703 (@20μM)

0.25 ul JDW 704 (@20μM)

0.125 ul Homemade TAQ

1.0 ul genomic DNA

PCR PROGRAM:

Step 1: 95°C – 3 minutes

Step 2: 94°C – 30 seconds


Step 3: 63°C – 45 seconds

Step 4: 72°C – 45 seconds

Step 5: GO TO STEP 2, x 35 cycles

Step 6: 72°C X 5 minutes

Step 7: 16°C FOREVER



*Run on a 2.5% gel

From the paper:

Primer sequences for Alk1 alleles are as follows: 3F cagcacctacatctgggtggaga; 3R actgttctctcggagcctgtc; 6F cctggacagcgactgtactac; and 6R gcccattgctctctcaaac. Combinations of 3F and 3R primers for the 3flox allele detection, 6F and 6R primers for the 2flox allele detection, and 3F and 6R primers for the 1flox (null) allele detection.